

Allergen, White Pine (t16) IgE

Test ID: 3706

CPT: 86003

Clinical Significance:

Latin name: Pinus strobus

Source material: Pollen

Family: Pinaceae

Common names: White pine, Eastern white pine, Northern white pine, Weymouth pine

Not to be confused with the Australian pine tree (Casuarina equisetifolia) t73.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: