Allergen, Banana, (f92) IgE

Test ID:

9894

CPT:

86003

Clinical Significance:

A food, which may result in allergy symptoms in sensitised individuals.

Allergen Exposure

Banana is the common name for a fruit, and also for the herbaceous plants of the genus Musa that produce this popular fruit. Banana plants are of the family Musaceae. They are cultivated primarily for their fruit, and to a lesser extent for the production of fibre and as ornamental plants.



Container:	
Serum Separator Tube (SST [®])	

Transport Temperature:

Room Temperature

Specimen:	
Serum	

Specimen Stability:
Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:
Mon-Sat

Collection Instructions: