

AA- Allergy, Seeds Panel

Test ID: 99984

CPT: 82785 , 86003 X 5

Clinical Significance:

Profile Components:

Sesame, f10, IgE
Poppy seed, f224, IgE
Flaxseed, rf333, IgE
Cotton seed, K83, IgE
Sunflower Seed, k84, IgE



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: