

# Allergen, Beef, (f27) IgE

Test ID: 3030

CPT: 86003

## Clinical Significance:

Several allergens present in cow dander and hair are also found in beef and milk.

Allergen Exposure:

Meat from bovines.

Unexpected exposure:

Undeclared ingredient in processed meat, like sausages, patês, etc.

Bacterial vaccines containing bovine serum albumin.

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: