

Allergen, Beta-lactoglob, (f77) IgE

Test ID: 3518

CPT: 86003

Clinical Significance:

Source material: nBos d 5 is purified from milk extract .
Milk contains more than 40 proteins, and all of them may act as human species antigens.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: