

Allergen, Birch,(t3) IgE

Test ID: 3524

CPT: 86003

Clinical Significance:

Source material: Pollen

Common names: Common Silver Birch, Common Birch, Birch, Birch tree

Synonym: *B. pendula*

Birches, *Betula* spp., are boreal and temperate trees and shrubs, with around 50 species worldwide and 12 species native to North America (1). River Birch (*B. nigra*) and White or Paper Birch (*B. papyrifera*) are common to North America, whereas Common Silver Birch (*B. verrucosa*) is the prevalent European species (2).

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: