## Allergen, Blueberry (f288) IgE

•	D:	I	st	Ге	٦
---	----	---	----	----	---

### 2568

CPT:

86003

#### **Clinical Significance:**

Source material: Frozen fruit

Common names: Blueberry, European blueberry, Lowbush blueberry, Highbush blueberry, Whinberry, Whortleberry, Bilberry

#### **Profile Components:**

# Labtech Diagnostics

	Container:
ıbe (SST®)	Serum Separa

Transport Temperature:

Room Temperature

Specimen Stability:
Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:
Mon-Sat

**Collection Instructions:**