

Allergen, Brazil Nut, (f18) IgE

Test ID: 3084

CPT: 86003

Clinical Significance:

Latin name: Bertholletia excelsa

Source material: Shelled nuts

Family: Lecythidaceae

Common names: Brazil nut, Para-nut, Cream nut

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: