
## Allergen, Buckwheat (f11) IgE

Test ID:	2811	CPT:	86003
Clinical Significance:			
Latin name: Fagopyrum e	sculentum		
Source material: Whole se	eed		
Family: Polygonaceae			
Common names: Buckwh	eat, Beechwheat, Fagopyrum, French	wheat, Garden buckwheat	



Conta	iner:		
Serum	Separator 1	ube (SST®)	
	•	, ,	

Transport Temperature:

Room Temperature

Specimen:			
Serum			

Specimen Stability:		
Room temperature: 28 days		
Refrigerated: 28 days		
Frozen: 1 year		

Reject Criteria:

<b>Days Perfor</b>	med:
Mon-Sat	

**Collection Instructions:**