

Allergen, *C. albicans*, (m5) IgE

Test ID: 3088

CPT: 86003

Clinical Significance:

Latin name: *Candida albicans*
Source material: Disrupted cells
Yeast

Candida-specific IgE antibodies have been demonstrated in asthma and rhinitis (1).

C. albicans is common in soil, organic debris and in humans where it occurs as a saprophyte in the nasopharynx and feces.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: