

Allergen, Carrot (f31) IgE

Test ID: 2831

CPT: 86003

Clinical Significance:

Latin name: Daucus carota
Source material: Fresh frozen juice
Family: Apiaceae
Common names: Carrot

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: