Allergen, Casein, (f78) IgE

lest ID:	3516	CPT:	86003
Clinical Significance:			
Latin name: Bos spp.			
Source material: Purified be	ovine casein		
Casein is a major allergen in	n milk and the main protein consti	tuent in cheese.	
Allergen Exposure			
Milk and dairy products, es	pecially cheese.		
Food containing milk.			
Casein is a heat stable majo	or allergen in milk and the main pro	otein constituent in cheese.	

Profile Components: Labtech Diagnostics

Container:	Transp
Serum Separator Tube (SST [®])	Room

Transport Temperature: Room Temperature

Specimen:	
Serum	

Specimen Stability:		
Room temperature: 28 days		
Refrigerated: 28 days		
Frozen: 1 year		

Reject Criteria:

Days Performed:	
Mon-Sat	

Collection Instructions: