

Allergen, Chicken, (f83) IgE

Test ID:	3506	CPT:	86003
Clinical Significance:			
Latin name: Gallus spp.			
Source material: Freeze-dri	ed, raw meat		
Family: Galliformes			
Chicken serum albumin ma	y induce a cross-sensitization betwe	een egg and chicken proteins.	





Transport Temperature:

Room Temperature

Specimen:			
Serum			

Specimen Stability:		
Room temperature: 28 days		
Refrigerated: 28 days		
Frozen: 1 year		

Reject Criteria:

Days Performed:	
Mon-Sat	

Collection Instructions: