Allergen, Clam, F207, IgE

9944 Test ID: CPT: 86003 **Clinical Significance:** Latin name: Clam Source material: Fresh frozen muscle Family: Veneridae Clam is closely related to quahogs and cockles. Allergen Exposure Cooked, often in sauces. Clinical Experience Occupational asthma and food allergy caused by clam have been reported **Profile Components:** Labtech Diagnostics **Container:** Transport Temperature: Specimen: Serum Separator Tube (SST®) Room Temperature Serum Specimen Stability: Reject Criteria: **Days Performed:** Room temperature: 28 days Mon-Sat Refrigerated: 28 days Frozen: 1 year **Collection Instructions:**