

Allergen, Cocklebur (w13) IgE

Test ID: 2413

CPT: 86003

Clinical Significance:

Latin name: Xanthium commune

Source material: Pollen

Family: Asteraceae

Common names: Cocklebur, Rough Cocklebur, Common Cocklebur

Allergen Exposure

Cocklebur is native to Europe, Asia, southern Canada, and most of the United States, Mexico and Central America. Cocklebur is now found worldwide. Another species that is less common but widespread across North America is called Spiny Cocklebur

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: