Allergen, Coconut, (f36) IgE

CDT

2024

163(10).	3034	CF1. 0000
Clinical Significance:		
Latin name: Cocos nucifera		
Source material: Fresh coconut meat		
Family: Arecaceae		
Common names: Coconut, Common Coc	conut	
There have in the past been 60 other spe	ecies classified under the genus Cocos, bu	t the Coconut palm is actually monotypic,
meaning that within the genus Cocos on	ly one species, nucifera, is now recognise	d, though within this species are numerous

Profile Components: Labtech Diagnostics

Container:	
Serum Separator Tube (SST [®])	

Tost ID.

varieties.

Transport Temperature:

Room Temperature

Specimen:		
Serum		

86003

Specimen Stability:
Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:
Mon-Sat

Collection Instructions: