
## Allergen, Common Pigweed, (w14) IgE

lest ID:	3536	CPT:	8600
Clinical Significance:			
Latin name: Amaranthus retr	oflexus		
Source material: Pollen			

Family: Amaranthaceae

Common names: Common Pigweed, Redroot Pigweed

Note: Lamb's quarters (Chenopodium album) is occasionally also called Pigweed or Smooth pigweed but does not belong to the Amaranthaceae family. There is a particular resemblance in the cotyledon stage, but Lamb's quarters cotyledons often have a mealy grey cast and the first true leaves are alternate, unlike those of any of the Pigweed species.



Contair	ner:			
Serum	Separato	r Tube (	(SST®)	

Transport Temperature:

Room Temperature

Specimen Stability:	
Room temperature: 28 days	
Refrigerated: 28 days	
Frozen: 1 year	

Reject Criteria:

Days Performed:
Mon-Sat

**Collection Instructions:**