
Allergen, Corn, (f8) IgE

Test ID:	3006	CPT:	86003
Clinical Significance:			
Code: f8			
Latin name: Zea mays			
Source material: Untreat	ed planting seeds		
Family: Poaceae (Gramin	eae)		
Common names: Maize,	Corn, Sweet Corn, Indian Corn, Field C	Corn	
Cultivars within the genu	s may be divided into 6 general types	: Popcorn (everta), Flint corn (indurata), De	ent corn (indenta), Flour
corn (amylacea), Sweet c	orn (saccharata) and Pod corn (tunica	ata).	

There are, however, only 2 basic types: "Sweet corn" is distinguished from "Field corn" by the high sugar content of the kernels



Container:	Transport Temperature:	Specimen:
Serum Separator Tube (SST [®])	Room Temperature	Serum
· · · · ·	·	

specifien Stability.		
Room temperature: 28 days		
Refrigerated: 28 days		
Frozen: 1 year		
Refrigerated: 28 days		

Reject Criteria:

Days Performed:
Mon-Sat

Collection Instructions: