Allergen, Cottonwood (t14) IgE

Test ID: 3204 CPT: 86003 **Clinical Significance:** Latin name: Populus deltoides Source material: Pollen Family: Salicaceae Common names: Cottonwood, Poplar tree The botanical family willow (Salicaceae) consists of Populus and Salix, the latter being true willows. Populus deltoides is the common designation for species of trees in this genus, but in fact Populus can be divided into 3 distinct groups: cottonwoods, aspens and poplars. **Profile Components:** Labtech Diagnostics **Container:** Transport Temperature: Specimen: Serum Separator Tube (SST®) Room Temperature Serum Specimen Stability: Reject Criteria: **Days Performed:** Room temperature: 28 days Mon-Sat Refrigerated: 28 days Frozen: 1 year

Collection Instructions: