

Allergen, Cranberry (rf341) IgE

Test ID: 38254

CPT: 86003

Clinical Significance:

Latin name: *Vaccinium oxycoccus*

Source material: Fresh fruit

Family: Ericaceae

Common names: Bog cranberry, Small cranberry, Bigger cranberry, Wild cranberry, Swamp cranberry, Marshwort, Fenne berry, Marsh whortleberry, Bounceberry, Craneberry

Synonyms: *V. microcarpos*, *V. palustre*, *V. hagerupii*, *Oxycoccus quadripetala*, *O. palustris*, *O. hagerupii*, *O. intermedius*, *O. microcarpus*, *O. ovalifolius*, *O. oxycoccus*, *O. quadripetalus*

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: