

# Allergen, *Curvularia lunata* (m16) IgE

Test ID: 6680

CPT: 86003

## Clinical Significance:

Latin name: *Curvularia lunata*

Source material: Spores and mycelium

There appears to be extensive cross-reactivity between *Stemphylium* and *Curvularia* and *Alternata* (3).

Allergen Exposure

*Curvularia* is a facultative pathogen and may cause leaf spots and seedling blight. It is also seen on castor beans, cotton, rice, barley, wheat, and corn.

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: