

# Allergen, *D farinae*, (d2) IgE

Test ID: 3718

CPT: 86003

## Clinical Significance:

Code: d2  
Latin name: Dermatophagoides farinae  
Source material: Whole body culture  
Family: Pyroglyphidae  
Common names: House dust mite, dust mite

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days  
Refrigerated: 28 days  
Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: