

Allergen, Dill (f277) IgE

Test ID: 2918

CPT: 86003

Clinical Significance:

Latin name: Anethum graveolens

Family: Apiaceae

Common names: Dill, False Anise, Bastard Fennel, Russian Parsley, Swedish Parsley

Synonym: Peucedanum graveolens

Spice

A spice, which may result in allergy symptoms in sensitised individuals.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: