## Allergen, Alder (t2) IgE

Test ID.		Tes	t I	D:
----------	--	-----	-----	----

3735


CPT:

86003

## **Clinical Significance:**

The Fagales order consists of trees from the Betulaceae family (grey alder and birch tree), along with hazel, hornbeam, oak, and chestnut tree.

Grey alder is the most widely distributed alder in Europe and western North America.



Contai	ner:			
Serum	Separate	or Tube	(SST®)	

Transport Temperature:

Room Temperature

Specimen:		
Serum		

Specimen Stability:
Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed: Mon-Sat		
Mon-Sat		

**Collection Instructions:**