

Allergen, Duck Feathers (e86) IgE

Test ID: 2664

CPT: 86003

Clinical Significance:

Latin name: *Anas platyrhynca*

Source material: Feathers

Family: Anatidae

Common names: Duck; also see below

Direct or indirect contact with bird allergens may cause sensitisation. Bird allergens may be major components of house dust.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: