

Allergen, Fire Ant, (i70) IgE

Test ID:	3544	CPT:	86003
Clinical Significance:			
Latin name: Solenopsis	invicta		
Source material: Whole	insect bodies		
Family: Myrmicidae			
Insect			
An insect, which may re	sult in allergy symptoms in sensitised	individuals.	





Room Temperature

Transport Temperature:

Specimen Stability:		
Room temperature: 28 days		
Refrigerated: 28 days		
Frozen: 1 year		

Reject Criteria:

Days	Performed:
N 4	C - +

Mon-Sat

Collection Instructions: