Allergen, Firebush (w17) IgE

Test ID: 2417 CPT: 86003

Clinical Significance:

Latin name: Kochia scoparia (Synonyms: Bassia scoparia, Chenopodium scoparia)

Source material: Pollen

Family: Amaranthaceae (Chenopodiaceae)

Common names: Firebush, Kochia, Common Kochia

Allergen Exposure

Kochia is native to southern and eastern Russia, Europe and Asia. It is now naturalised across the northern half of the United States and spreading south-westwards. It is found in many other areas of the world. Kochia is a major source of pollen.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: