

Allergen, *Fusarium moniliforme* (m9) IgE

Test ID: 3729

CPT: 86003

Clinical Significance:

Latin name: *Fusarium proliferatum*/F.moniliforme
Source material: Spores and mycelium
Family: Hypocreaceae
Common names: *Fusarium proliferatum*, *Fusarium moniliforme*
Mould
A mould, which may result in allergy symptoms in sensitised individuals.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: