

Allergen, Goat Epithelia (e80) IgE

Test ID: 2656

CPT: 86003

Clinical Significance:


Latin name: Capra hircus

Family: Bovidae

Common names: Goat

Direct or indirect contact with animal allergens frequently causes sensitisation.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: