

Allergen, Goose Feathers (e70) IgE

Test ID: 2661

CPT: 86003

Clinical Significance:

Latin name: Anser anser

Source material: Feathers and serum proteins

Family: Anatidae


Common names: Goose, Toulouse Goose, Grey Goose, Greylag

The genera Anser and Branta should both be considered. Anser includes many wild Geese with well-known domestic relatives.

This group includes Anser anser (European) and species of Africa and Asia. Branta includes wild Geese like the Canada Goose (Branta Canadensis) and the Brants (most notably Branta bernicla).

Direct or indirect contact with bird allergens may cause sensitisation. Bird allergens may be major components of house dust.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: