

Allergen, Grape (f259) IgE


Test ID: 2675

CPT: 86003

Clinical Significance:

Latin name: Vitis vinifera
Source material: Fresh fruit
Family: Vitaceae
Common names: Grape
Food
A food, which may result in allergy symptoms in sensitised individuals.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: