

Allergen, Green Bean, (f315) IgE

Test ID: 3554

CPT: 86003

Clinical Significance:

Latin name: Phaseolus vulgaris

Source material: Fresh beans

Family: Fabaceae (Leguminosae)

Common names: Green bean, Common bean, French bean, String bean, Snap bean, Wax bean, Haricot bean

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: