

# Allergen, Guinea Pig Epithelia (e6) IgE

Test ID: 2606

CPT: 86003

## Clinical Significance:

Latin name: *Cavia porcellus*

Source material: Epithelium and dander

Family: Caviidae

Common names: Guinea pig, Cavy

Direct or indirect contact with animal allergens frequently causes sensitisation.

Animal allergens are major components of house dust and laboratory dust.

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: