Allergen, Hazelnut, (f17) IgE

3080 Test ID: CPT: 86003 **Clinical Significance:** Latin name: Corylus avellana Source material: Shelled nuts Family: Betulaceae (Corylaceae) Common names: Hazel nut, Hazelnut, Filbert, Cobnut, Cob The terms "Filbert" and "Hazel nut" are often used interchangeably for nuts from all plants in the genus Corylus, such as C. silvestris, C. maxima and C. colurna. **Profile Components:** Labtech Diagnostics **Container:** Transport Temperature: Specimen: Serum Separator Tube (SST®) Room Temperature Serum Specimen Stability: Reject Criteria: **Days Performed:** Room temperature: 28 days Mon-Sat Refrigerated: 28 days Frozen: 1 year **Collection Instructions:**