
## Allergen, Herring (f205) IgE

Test ID:	2707	СРТ:	86003
Clinical Significance:			
Latin name: Clupea hare	ngus		
Source material: Fish mu	iscle		
Family: Clupeidae			
Small herrings are some	times canned and sold as sardines.		
Allergen Exposure			
Herring is eaten raw, co	oked, pickled, fried or smoked. Canned	"sardines" may sometimes be herrings.	



Conta	iner:			
Serun	n Separato	r Tube (S	SST®)	
	•	•		

Transport Temperature:

Room Temperature

Specimen:		
Serum		

Specimen Stability:		
Room temperature: 28 days		
Refrigerated: 28 days		
Frozen: 1 year		

Reject Criteria:

Days Performed:
Mon-Sat

**Collection Instructions:**