

Allergen, Honey bee, (i1) IgE

Test ID: 3300

CPT: 86003

Clinical Significance:

Latin name: Apis mellifera


Source material: Venom

Family: Apidae

Order: Hymenoptera

Honeybee is the only stinging Hymenoptera that nearly always leaves its fluked stinger in the skin of the victim.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: