

Allergen, Horse Dander (e3) IgE

Test ID: 2603

CPT: 86003

Clinical Significance:

Latin name: Equus caballus

Source material: Dander

Family: Equidae

Common names: Horse, Domestic Horse

Direct or indirect contact with animal allergens frequently causes sensitisation.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: