

AA- Allergy, Meat Panel

Test ID: 99987


CPT: 82785, 86003 x 5

Clinical Significance:

See individual tests

Profile Components:

Turkey Meat, f284, IgE
Lamb, f88, IgE
Beef, f27, IgE
Pork, f26, IgE
Chicken, F83, IgE



Labtech Diagnostics

Container:

2 X Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: