

# ANA DIRECT SCREEN, MULTIPLEX IMMUNOASSAY

Test ID: 250

CPT: 86038

## Clinical Significance:

Detect antibodies to nuclear antigens using Multiplex flow immunoassay

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 4 days  
Refrigerated: 7 days  
Frozen: 30 days

### Reject Criteria:

Gross hemolysis • Grossly lipemic  
• Grossly icteric

### Days Performed:

Mon-Fri

## Collection Instructions: