

Aspergillus Antibodies, Immunodiffusion

Test ID: 20341

CPT: 86606 (x3)

Clinical Significance:

Profile Components:

Aspergillus flavus, Aspergillus fumigatus, Aspergillus niger



Labtech Diagnostics

Container:

Serum Separator Tube (SST)

Transport Temperature:

Refrigerated

Specimen:

1 mL Serum

Specimen Stability:

Room temperature: 2 days
Refrigerated: 14 days
Frozen: 1 year

Reject Criteria:

Body Fluid Samples
Avoid Repeated Freeze/Thaw Cycles
Contaminated
Hemolysis
Lipemia

Days Performed:

Thu

Collection Instructions:

1 mL serum. Allow SST to clot in an upright position for at least 30 minutes, then centrifuge sample within 2 hours of collection. Refrigerate