

# Anti-Centromere B Antibodies


Test ID: 8220

CPT: 86038

## Clinical Significance:

Centromere B Antibody is diagnostic for the form of scleroderma known as CREST (calcinosis, Raynaud's phenomenon, esophageal immotility, sclerodactyly, and telangiectasia).

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 4 days  
Refrigerated: 7 days  
Frozen: 30 days

### Reject Criteria:

Gross hemolysis • Gross lipemia •  
Grossly icteric

### Days Performed:

Mon-Fri

## Collection Instructions: