

Allergy, Weeds Panel

Test ID: 99997

CPT: 86003

Clinical Significance:

See individual tests

Profile Components:

Giant Ragweed, w3, IgE
Mugwort, w6, IgE
English plantain, w9, IgE
Lamb's Quarters (Goose Foot) ,w10, IgE
Common Ragweed, w1, IgE
Pigweed Rough, w14, IgE
Nettle,w20, IgE
Sheep Sorrel, w18, IgE

Rough Marsh Elder, w16, IgE
Firebush, w17 IgE
Cocklebur, w13, IgE
Goldenrod, w12, IgE



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)X2

Transport Temperature:

Room temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: