

# Amphetamine, LCMS Confirmation, Ur

Test ID: 850823

CPT: G0480

## Clinical Significance:


Toxicology LCMS confirmation are used to identify substances affecting a patient and to guide the clinician to predict future toxic effects, to confirm the differential diagnosis or to guide therapy. Accurate diagnosis of clinical intoxication secondary to illicit drug use based on clinical history and physical examination may be difficult without laboratory confirmation. This testing is also necessary when multiple drug ingestion is involved, as the effects of one drug may mask the clinical signs and symptoms of the effects of other drugs.

## Profile Components:

Amphetamine

Present in the following medication :

Adderall, Adzenys XR-ODT, Benzedrine, Dexedrine, Dyanavel XR, Evekeo, Vyvanse, ProCentra, Zenzedi



# Labtech Diagnostics

### Container:

Sealed Urine Container

### Transport Temperature:

Room temperature

### Specimen:

Urine

### Specimen Stability:

Room temperature: 10 days

Refrigerated: 30 days

Frozen: 30 days

### Reject Criteria:

### Days Performed:

Mon-Fri

## Collection Instructions:

Urine temperature monitoring cup is recommended for samples integrity .