

Alpha-Fetoprotein (AFP) and AFP-L3

Test ID: 19529

CPT: 82107

Clinical Significance:

This assay is intended for use in the assessment of risk for the development of hepatocellular carcinoma (HCC) in patients with chronic liver disease.

Profile Components:



Labtech Diagnostics

At the Forefront of
Clinical Laboratory Medicine

Container:

Serum Separator Tube (SST)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 72 hours
Refrigerated: 5 days
Frozen: 2 years

Reject Criteria:

Hemolysis • Lipemia

Days Performed:

Days Performed: Thursday AM shift
Turn Around Time: 1-8 Business Days

Collection Instructions:

Separate serum from cells within 45 minutes of collection.