

# Alcohol, EtG, Qualitative, Random Urine

Test ID: 90418

CPT: 80321

## Clinical Significance:

A positive urine alcohol is only indicative of recent use and cannot be used to determine impairment. This test should be restricted to monitoring patients in drug treatment programs where any alcohol use is prohibited

## Profile Components:



# Labtech Diagnostics

### Container:

10 mL random urine

### Transport Temperature:

Room temperature

### Specimen:

Urine

### Specimen Stability:

Room temperature: 14 days  
Refrigerated: 14 days  
Frozen: 30 days

### Reject Criteria:

Urine with preservative

### Days Performed:

Mon-Sat

## Collection Instructions:

Random urine collection