

Allergy Environmental Panel,IgE

Test ID: 99886

CPT: 86003 X 16

Clinical Significance:

Profile Components:

Dog Dander, E5, Cat Dander, E1, Lamb's Qtrs, W10, Common Ragweed, W1, Bermuda grass, G2, Johnson Grass, G10, June Grass, KYB, G8, Cockroach, I6, Alternaria tenuis, M6, Aspergillus fumigatus, M3, Elm American, T8, Oak, T7, Pecan Tree, T22, Maple, Box Elde, T1, D. Pteronyssinus, D1, C. herbarum, M2.



Labtech Diagnostics

Container:

Serum Separator
Tube (SST)

Transport Temperature:

Room temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Collection Instructions:

Days Performed:

Mon-Sat

Reject Criteria