

# A. PHAGOCYTOPHILUM AB, IgM

Test ID: C2778

CPT: 86666

## Clinical Significance:

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST)

### Transport Temperature:

Refrigerated

### Specimen:

Serum

### Specimen Stability:

after separation from cells  
Room Temperature: 48 Hours  
Refrigerated: 2 Weeks  
Frozen: 1 Year

### Reject Criteria:

hemolysis  
lipemia  
icterus  
samples with bacterial contamination  
turbid samples  
heat-inactivated samples

### Days Performed:

## Collection Instructions:

1 mL serum. Allow SST to clot in an upright position for at least 30 minutes, then centrifuge sample within 2 hours of collection. Refrigerate. Note: parallel testing is preferred and convalescent sample must be received within 30 days of receipt of acute sample. Clearly label tube as acute or convalescent.