

Apolipoprotein B/A Ratio

Test ID: 4109

CPT: 82172 X2

Clinical Significance:

Evaluation of patients with signs and symptoms compatible with connective tissue diseases
Parietal Cell Antibody

Profile Components:



Labtech Diagnostics

Container:

Serum Specimen Tube (SST)

Transport Temperature:

Refrigerated

Specimen:

Serum

Specimen Stability:

Room temperature: 7 days
Refrigerated: 10 days
Frozen: 90 days

Reject Criteria:

Avoid Repeated Freeze/
Thaw Cycles

Days Performed:

Mon-Sat

Collection Instructions:

1 mL serum. Allow SST to clot in an upright position for at least 30 minutes, then centrifuge sample within 2 hours of collection. Refrigerate.